

FORA

FITNESS SCALE

W310

Dimension (inches) : 13.2 x 13.2 x 1.44"
Weight capacity: 180kg / 396 lb / 28 st:13lb
Auto-on/off function
Resolution : 0.1kg / 0.2lb / 0.2lb
Weighing Unit: kg/ lb/ st



- Measure your weight and BMI
- Large LCD display
- Built-in 10 user profiles
- Transmit your data from the W300 to the FORA TeleHealth system for even greater data management.

FORA[®]

FITNESS SCALE

W320

Dimension (inches) : 13.2 x 13.2 x 1.44"
Weight capacity: 180kg / 396 lb / 28 st:13lb
Auto-on/off function
Resolution : 0.1kg / 0.2lb / 0.2lb
Weighing Unit: kg/ lb/ st



- Large LCD display
- Large memory
- Transmit your data from the W320 to the FORA TeleHealth system for even greater data management