

# Oscilla® SM950 Diagnostic Memory Audiometer

---



# Oscilla® SM950 Diagnostic Memory Audiometer

## Features:

- 11 frequencies: 125 - 8000 Hz (selectable)
- -10 dB to 110 dB in 5 dB
- Narrow band masking
- USB port for PC-connection

## AudioConsole software\*:

- User-friendly software solution for on-screen audiometry
- Fast and easy PC-connection via USB (plug & play)
- Settings can be customised according to user preferences
- Unlimited database and storage capacity
- Generate PDF reports
- Noah compatible

\* AudioConsole is an optional accessory for SM950

## Manual and automatic tests:

- Manual test
- 20 dB automatic test
- 20 dB RANDOM automatic test
- The Hughson Westlake automatic test
- Weber test (only in Oscilla® AudioConsole)
- Ear protection test (only in Oscilla® AudioConsole)

## Portable and user-friendly:

- Simple operation
- Portable and compact in size
- Noise-free push buttons

**Standard accessories:** SM950 audiometer, patient headphones, patient response, bone conductor, power supply

**Optional extras:** AudioConsole Software, USB cable, carrying bag

**Dimensions:** 295 x 180 x 55 mm **Power supply:** AC 100-240V, 50-60Hz to 12V DC  
**Weight:** 1000 g (device only)

Features and specifications are subject to change without prior notice

